

Latest advancements in Endodontics- I



Dr Saaid Al Shehadat

Year 4
8/10/2019

Endodontics Progress at 4 levels

- Tools : diagnosis, RCT procedures
- Materials: irrigation, obturation and restoration
- Techniques
- Concepts




Tools



Diagnosis

- Pulp vitality test
- CBCT



Pulp tests

Sensibility tests

- Cold Test
- Heat Test
- Electrical Test
- Other Tests: cavity test, anesthetic test

Shortcomings of sensibility test

- Accompanied by un-pleasant sensations for the patient
- false-positive response, the nervous system may remain reactive, even if all the surrounding tissues are degenerated
- false-negative response, in cases of carious (degenerescent), teeth with recent traumatism, or teeth with open apices

Recently available pulp vitality tests

- Laser Doppler flowmetry (LDF)
- Pulp oximetry
- Measurement of temperature of tooth surface
- Transillumination with fiber-optic light
- Dual wavelength spectrophotometry
- Plethysmography
- Detection of interleukin—1 beta
- Xenon—133
- Hughes probeye camera
- Gas desaturation
- Radiolabeled microspheres
- Electromagneticflowmetry


Cohen 2016

Tools

★ Diagnosis

Cone Beam Computed Tomography (CBCT)

Cone Beam Volumetric Tomography (CBVT)




CBCT

- Uses a cone-shaped beam of radiation to acquire a volume in a single 360-degree rotation
- The volume acquired by a CBCT is composed of voxels
- A **voxel** is a 3-D **pixel**



CBCT

With the help of viewer software, the clinician is able to scroll through the entire volume and simultaneously view **axial, coronal, and sagittal** 2-D sections that range from **0.125–2.0 mm** thick



CBCT systems

limited (dental or regional)

- The scan volume: 40 –100 μm
- Voxel: 0.1–0.2 mm^3

full (ortho or facial)

- The scan volume: 100 –200 μm
- Voxel: 0.3–0.4 mm^3

limited CBCT systems offer higher resolution and are better

Cotton, T.P., Geisler, T.M., Holden, D.T., Schwartz, S.A. and Schindler, W.G., 2007. Endodontic applications of cone-beam volumetric tomography. *Journal of endodontics*, 33(9), pp.1121-1132.

CBCT systems

Limited CBCT is preferred for most **endodontic** applications:

1. Increased resolution, for (calcified/accessory canals, missed canals)
2. Focus on anatomical area of interest
3. Decreased radiation exposure to the patient
4. Time savings due to smaller volume to be interpreted
5. Smaller area of responsibility



CBCT


CBCT has become increasingly important in treatment planning and diagnosis for

- Implant
- Oral surgery
- Orthodontics
- Endodontics

CBCT in Endodontics

A drawback of periapical X-ray

The interpretation can be confounded due to regional anatomy as well as **superimposition** of both the teeth and surrounding structures




- Curvature (buccal/ lingual)
- Exact length
- Extension of lesions

Potential endodontic applications

- diagnosis of endodontic pathosis
- assessment of pathosis **of non-endodontic origin**
- calcified/accessory canals, missed canals
- canal morphology
- evaluation of root **fractures** and trauma
- analysis of external and internal root **resorption**
- **invasive cervical resorption**
- pre-surgical planning
- to determine the nature of the peri-radicular lesion (granuloma vs cyst)

CBCT in Endodontics

Missing canal




Cotton et al., 2007

CBCT in Endodontics

Nasopalatine duct cyst


- RCT of 11
- internal resorption on #11?
- 2 PA: radiolucency shifts entirely off the root structure
- external resorptive defect? pathosis within the palatal bone?
- CBCT:



Cotton et al., 2007

CBCT in Endodontics


Internal resorption



Cotton et al., 2007


CBCT in Endodontics

Pre-surgical anatomic assessment




Cotton et al., 2007

CBCT in Endodontics



Pre-operative

Calcium hydroxide dressing



Post-operative




6 months re-call

<https://endolit.com/endoappforum/>

CBCT in Endodontics

Pre-operative

6 months recall




<https://endolit.com/endoappforum/>

CBCT in Endodontics

Dynamic Guided Endodontic

Digital impressions + CBCT scans + 3D printing

To negotiate calcified canal without perforation



van der Meer et al. (2016)

Tools


★ Root Canal Preparation

- Access cavity
- Irrigation
- Root canal shaping



Access cavity and canal instrumentation

Round burs → Conical carbide bur



A new concept: Minimal Invasive Endodontics

Access cavity and canal instrumentation


Conical carbide bur

- Self centering
- Safer
- Minimal invasive




Access cavity and canal instrumentation

round burs



conical carbide bur




Access cavity and canal instrumentation

Ultrasonic tips (US)

- access refinement
- finding calcified canals
- removal of attached pulp stones
- removal of posts
- removal of broken instruments



Access cavity and canal instrumentation



Access cavity and canal instrumentation

Terauchi file removal kit (TFRK)



Tools


★ Root Canal Irrigation



Modern concepts in root canal therapy procedures

The complexity of the root canal

- Isthmus may connect the canals
- only 45-55% of canal walls are actually touched
- The manual irrigation is not sufficient
- The use of a high concentration of NaOCl or EDTA lead to only marginally better results



3D Micro CT: Canal morphology before instrumentation (green); canal walls touched using a standard NITI file (red).
*Courtesy of Dr. Frank Pappal (Switzerland)

Sodium hypochlorite (NaOCl)

- Concentrations: 0.5% and 6%
- Both lower and higher concentrations are equally efficient in reducing the number of bacteria in infected root canals
- The time needed to inhibit bacterial growth and tissue dissolving effect of NaOCl irrigant is related to its concentration



Sodium hypochlorite (NaOCl)

Temperature

Increasing the temperature of hypochlorite irrigant to 60°C, significantly increases its antimicrobial and tissue-dissolving effects



Heating devices for NaOCl syringe

Irrigation Devices In Endodontic

Manual


- 1- Plastic syringe with needles
- 2- Brushes
- 3- Manual Dynamic Agitation

Machine Assisted


- | | |
|--|---------------------------------|
| 1- Rotary brushes | 4- Ultrasonic |
| 2- Continuous irrigation during rotary instrumentation | 5- Pressure alternating devices |
| 3- Sonic | 6- Recent advance system |

Plastic syringe with needles

Several needle-tip design have been introduced to facilitate effectiveness and minimize safety risks



Beveled Monoject Closed end



Flexiglide needle

Plastic syringe with needles

Navitip


Navitip tip is a small, flexible cannula easily delivers endodontic irrigants and sealers



Plastic syringe with needles

Navitip


- Mostly rigid shanks, with the last 5 mm flexible to facilitate navigation through curved canals
- 29- 30 gauge, lengths 17, 21, 25 and 27 mm



Plastic syringe with needles

Navitip


- ✓ Safe, can easily be inserted to apical third



Recent needles

Navitip FX

31g
Double Sideport Irrigator



www.youtube.com/watch?v=3FVXN1sCKf8

2- Brushes

Endobrush

- ✓ used only as an adjuncts
- ✓ Nylon bristles set in twisted wire

Disadvantages

- ✗ Can't be used till working length
- ✗ Dislodgement of radiolucent bristles



3- Manual Dynamic Agitation

Gutta- percha points agitation

To facilitate the penetration of solutions in the canal, gently move a well-fitting gutta-percha master cone up and down in short 2- to 3-mm strokes within an instrumented canal



Basic principles of cleaning


To clean the table we don't need to cut its surface, but touch ALL table surface using some chemical and some force to remove dirt attached.



Root Canal Irrigation

Machine Assisted Devices

1. Rotary Brushes



Ruddle brush

Canal brush

2. Continuous Irrigation During Rotary Instrumentation

Quantec-E

✓ Self contained fluid delivery unit

✓ uses a pump console, 2 irrigation reservoirs & tubing



(Quantec-E)

3- Sonic devices

✓ combines battery-driven vibrations (9000- 10,000 cpm)



Rispisonic file



Vibringe


Endoactivator



4- Ultrasound

i- Active Ultrasonic Irrigation (AUI)

It is the simultaneous combination of ultrasonic irrigation and instrumentation. it has been almost discarded in the clinical



4- Ultrasound

ii- Passive Ultrasonic Irrigation (PUI)

- Frequencies of 25-30 kHz
- The mechanical energy warms the irrigant solution (NaOCl) and dislodges debris from canals




4- Ultrasound

ii- Passive Ultrasonic Irrigation (PUI)

- ✓ The energy is transmitted from an oscillating file to the irrigant in the root canal by means of ultrasonic waves
- ✓ The latter induces acoustic streaming and cavitation of the irrigant



Acoustic streaming



cavitation

4- Ultrasound

ii- Passive Ultrasonic Irrigation (PUI)

- Tips are small, parallel-shaped and non-cutting and can be used in the complete root canal
- 1-3 minutes of irrigation at the end of the root canal preparation (inserted 1mm short of working length)

Irri-Safe Tips

Clinical Research Dental Supplies




5- Pressure Alternation Devices

i- EndoVac System



5- Pressure Alternation Devices

i- EndoVac System




- Macro Cycle 30 seconds NaOCl (up and down)
- Micro Canula to WL: 30 seconds NaOCl
30 seconds EDTA
30 seconds NaOCl

<http://www.sybronendo.com/index/endovac-irrigation-technique>

5- Pressure Alternation Devices

i- EndoVac System

- Is based on a negative-pressure approach
- The irrigant placed in the pulp chamber is sucked down the root canal and back up again through a thin needle



(+) pressure
(traditional device)




(-) pressure (EndoVac)

5- Pressure Alternation Devices

i- EndoVac System

Advantages

- 1.Safety:** Less apical extrusion risk
- 2.Efficacy:** Better debridement 1 mm from working length
- 3.Success:** Negative apical pressure results in significantly less postoperative pain



5- Pressure Alternation Devices

i- EndoVac System

Relieve pressure from abscess



<http://www.sybronendo.com/index/sybronendo-clean-endovac-02>

5- Pressure Alternation Devices

ii- RinsEndo System

- ✓ Pressure phase: irrigant oscillating at a frequency of 1.6 Hz is drawn from an attached syringe and transported to the root canal via an adapted cannula.
- ✓ Suction phase (100 times per minute)

Disadvantage

- ✗ higher risk of apical extrusion of the irrigant




(Durr Dental Co)

5- Pressure Alternation Devices

ii- RinsEndo System

- Is based on a pressure-suction mechanism
- Contains:
 - * Titanium handpiece
 - * specially designed single-use cannulas
 - * A protective stopper guards against splatter and serves as a positioning device for the saliva ejector.



(Durr Dental Co)

6- Recent Advances in Irrigation System

i- LASERS

Nd:YAG Ho:YAG CO2 Er:YAG

- significantly reduced the number of bacteria
- It is inferior to NaOCl irrigation
- effective in removing and melting the smear layer



6- Recent Advances in Irrigation System

ii- Light-activated Disinfection (LAD)

- ✓ Photodynamic antimicrobial chemotherapy
- ✓ Photosensitizer (toluidine blue dye, methylene blue dye, etc)
- ✓ The canal is then filled with a photosensitizer and then illuminated with a light source (laser, white light, red light, or a light-emitting diode)



FotoSan: commercially available LAD

Thank You

